


### Wetness indicator

Each Abri-San diaper is provided with a wetness indicator, which clearly shows when the diaper should be changed. This also reduces the consumption of diapers.

### Hydrophobic

Standing cuffs of hydrophobic non-woven material offer, double protection against leakage.




### Plastic Backsheet

The foil is made of polyethylene which is converted into water and carbon dioxide when incinerated

### Colour codes

The diapers are provided with different backsheet colours to simplify identification.

## Security 24 hours a day



### TOP DRY®

The isolating Top Dry layer triples the distance between the skin and the part of the diaper where the liquid is absorbed. The diaper always feels dry and comfortable.

### Dry surface

The surface consists of a non-woven fabric that is permanent penetrable to liquids. So the diaper always feels pleasantly dry.

### Fast liquid dispersion - a dry feeling

The lengthwise channel system in the absorption core ensures a fast drain and dispersion of the liquid.

### Double absorption core

A super strong absorption core is placed in the middle of the diaper.

### Improved superabsorbent

The new superabsorbent not only keeps the liquid effectively in the diaper but also prevents the development of unpleasant odours.

